

Watch the
Total Concept!

Mobilyzer STANDARD

For the mobile patient:

- Internal medicine
- Orthopaedics
- Cardiology
- Surgery

Mobilyzer BASIC

For each patient who
needs support:

- Internal medicine
- Gynaecology
- Lung diseases
- Neurology
- Nephrology
- Rheumatology

Mobilyzer OPTION

For each patient who
needs support and who
carries equipment such
as infusion pumps:

- Oncology
- Cardiology
- Surgery



THE CONCEPT

The Mobilyzer concept consists of three variations:

1. **Mobilyzer Standard.**
2. **Mobilyzer Basic.**
3. **Mobilyzer Option.**

These variations can be customized by adding various accessories.

UNIQUE PRODUCT FEATURES

- A modular base on which accessories can be mounted.
- 5 wheels: Irregularities on the floor or thresholds can be passed easily because of the unique position of the fifth wheel.
- Smart storage: 21 units in an area of 2x2 meters.
- Ergonomics:
 - Ergonomic design for optimal use
 - Support through the handlebar
 - Height adjustable handlebar
- Non rotating handlebar by means of the oval shape of the upright tube of the frame.
- The Mobilyzer supports independence of the patient and encourages mobility.

SPECIFICATIONS

- Stainless steel frame (electropolished).
- Dimensions:
 - Height: max. 2000mm – 79.7"
 - Width: 600/650mm – 23.6"/25.6"
 - Length: 550/730mm – 21.7"/28.7"
- Weight capacity:
 - Max. handlebar: 130 kg – 287 lb
 - Max. infusion rod: 50 kg – 110 lb
 - Max. infusion hook: 25 kg – 55 lb
 - Max. capacity base: 250 kg – 551 lb
- Weight:
 - Mobilyzer Standard: ±10 kg- 22 lb
 - Mobilyzer Basic: ±12 kg- 26 lb
 - Mobilyzer Option: ±17,5 kg – 39 lb

