



**TechniCare**  
Improving Healthcare



# TOBA Shower chair

Shower chair for everyday use and complex care situations



## TOBA

The TOBA Showering chair is a stainless steel showering chair with an electrical adjustment, both in height as well in tilt. TechniCare delivers the TOBA complete, including head rest, foot rest and foldable arm rests. With this, the TOBA can be used directly.

## Complex Care situations

It is possible to add several accessories to make the TOBA suitable for complex care situations. For example TechniCare delivers custom made seatings, back rests and pads. The TOBA is often used both in hospitals, as well in home care.



# Specific aspects

- Electrical height and tilt adjustments
- Easy to maneuver in small bathrooms
- For every care situation a suited TOBA
- Stainless steel frame
- Easy to use remote control
- Tilt adjustment: -7° tot +32°
- Provided with emergency decrease
- Exchangeable battery



## Specifications TOBA

	basic	extra low	xl
Width Toba	60 cm	60 cm	72 cm
Width seating	42.5 cm	42.5 cm	55 cm
lowest position	48 cm	33 cm	48 cm
Highest position	106 cm	106 cm	106 cm
max. load	180 kg	180 kg	250 kg
Weight chair	50 kg	50 kg	50 kg
Tilt in lowest position	-7 tot 11	-7 tot 11	-7 tot 11
Tilt in highest position	16 tot 32	16 tot 32	16 tot 32



# Improving Healthcare

TechniCare's ambition is to improve care by providing technical innovation. TechniCare provides high quality products, developed in collaboration with the end user. We are a Dutch company with our own development and production facilities. For more information, visit [www.technicare.nl](http://www.technicare.nl)



Nikkelstraat 24  
2984 AM Ridderkerk

T +31 (0)85 0870 421  
[info@technicare.nl](mailto:info@technicare.nl)

[www.technicare.nl](http://www.technicare.nl)