



TechniCare
Improving Healthcare



Proceed Plaster Chair

Comfort and convenience



Compact plaster chair


Numerous orthopaedic technicians have been struggling with the same problem for years: large benches in small spaces that are often too small, and that don't offer a full range of functionalities. With the Proceed plaster chair/treatment chair, this is a thing of the past! The plaster chair is a compact chair with many adjustment possibilities that allow you to apply just about any plaster cast.

The plaster chair/treatment chair has been developed based on practical experience, which means that the chair has numerous practical features that contribute to enhanced user experience for both patient and orthopaedic technician.

Comfort and convenience

Comfort and convenience are the hallmark of the plaster chair - treatment chair, with patient comfort coming from the inclusion of not only nice cushions, but also comfortable entry and exit position. In addition, the plaster cast chair can be positioned so that any comfortable position for the patient can be easily adjusted.





Convenience for the orthopaedic technician is also imperative. The orthopaedic technician needs to be able to focus on healthcare, and the plaster chair should be of assistance to him/her in this regard. With intuitive operation and adjustment options, we ensure the easy application of a range of plaster casts.

Development

By and for orthopaedic technicians. The Proceed plaster chair has been developed in close cooperation with orthopaedic technician Leon Brouwers, who is also director of Medi-Chair. This chair is the result of many years of experience.

Accessoires



Knee support



Stainless steel material tray



Support cushion



Battery

Details

Width:	53 cm
Backrest adjustment:	85°
Height adjustment:	60-100 cm
Trendelenburg adjustment:	15°
Maximum weight:	150 kg
Wheels:	Ø100 mm
Material:	Steel coated

“Efficiency has never been so comfortable”

Mr Leon Brouwers, Owner of Medi-Chair and Orthopaedic Technician

Improving Healthcare

TechniCare aspires to improve healthcare by delivering technical innovation. TechniCare supplies high quality products, developed in collaboration with the end user. We are a Dutch company with our own development and production facilities. For more information see: www.technicare.nl



Roermondsestraat 54
5912 AK Venlo
The Netherlands

T +31 (0)6 403 94 738
info@medi-chair.nl

www.medi-chair.nl



Nikkelstraat 24
2984 AM Ridderkerk
The Netherlands

T +31 (0)85 0870 421
info@technicare.nl

www.technicare.nl